

Neurofilament heavy polypeptide

Catalog No: #AP86026



Package Size: #AP86026-1 50ug #AP86026-2 100ug #AP86026-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neurofilament heavy polypeptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	NF-H,200 kDa neurofilament protein,Neurofilament triplet H protein,Nfh
Accession No.	Uniprot ID:P16884
Uniprot	P16884
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nefh

Note: This product is for in vitro research use only