

Nuclear factor of activated T-cells, cytoplasmic 1

Catalog No: #AP86030



Package Size: #AP86030-1 50ug #AP86030-2 100ug #AP86030-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear factor of activated T-cells, cytoplasmic 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	NF-ATc1,NFAT transcription complex cytosolic component,NF-ATc
Accession No.	Uniprot ID:P98201
Uniprot	P98201
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NFATC1

Note: This product is for in vitro research use only