

## Nuclear factor NF-kappa-B p105 subunit

Catalog No: #AP86031



Package Size: #AP86031-1 50ug #AP86031-2 100ug #AP86031-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nuclear factor NF-kappa-B p105 subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DNA-binding factor KBF1,EBP-1,Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
Accession No.	Uniprot ID:Q63369
Uniprot	Q63369
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Nfkb1

Note: This product is for in vitro research use only