

Neurogenic locus notch homolog protein 2

Catalog No: #AP86033



Package Size: #AP86033-1 50ug #AP86033-2 100ug #AP86033-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurogenic locus notch homolog protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Notch 2, Motch B
Accession No.	Uniprot ID:O35516
Uniprot	O35516
GeneID	18129;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Notch2

Note: This product is for in vitro research use only