

Nuclear receptor subfamily 1 group I member 3

Catalog No: #AP86036



Package Size: #AP86036-1 50ug #AP86036-2 100ug #AP86036-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor subfamily 1 group I member 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Constitutive androstane receptor,CAR,Orphan nuclear receptor MB67,Car
Accession No.	Uniprot ID:O35627
Uniprot	O35627
GeneID	12355;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nr1i3

Note: This product is for in vitro research use only