

Nuclear receptor subfamily 5 group A member 2

Catalog No: #AP86037



Package Size: #AP86037-1 50ug #AP86037-2 100ug #AP86037-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor subfamily 5 group A member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	FTF/LRH-1,OR2.0
Accession No.	Uniprot ID:O42101
Uniprot	O42101
GeneID	395961;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NR5A2

Note: This product is for in vitro research use only