

NAD(P)H dehydrogenase [quinone] 1

Catalog No: #AP86042



Package Size: #AP86042-1 50ug #AP86042-2 100ug #AP86042-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD(P)H dehydrogenase [quinone] 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Azoreductase,DT-diaphorase,DTD,Menadione reductase,NAD(P)H:quinone oxidoreductase 1,Phylloquinone reductase,Quinone reductase 1,QR1,DIA4,NMOR1
Accession No.	Uniprot ID:P15559
Uniprot	P15559
GeneID	1728;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NQO1

Note: This product is for in vitro research use only