

Platelet endothelial cell adhesion molecule

Catalog No: #AP86047



Package Size: #AP86047-1 50ug #AP86047-2 100ug #AP86047-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet endothelial cell adhesion molecule
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PECAM-1,EndoCAM,GPIIA',PECA1,CD31
Accession No.	Uniprot ID:P16284
Uniprot	P16284
GeneID	5175;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PECAM1

Note: This product is for in vitro research use only