

Xaa-Pro dipeptidase

Catalog No: #AP86048



Package Size: #AP86048-1 50ug #AP86048-2 100ug #AP86048-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Xaa-Pro dipeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	X-Pro dipeptidase,Imidodipeptidase,Peptidase D,Proline dipeptidase,Prolidase,PRD
Accession No.	Uniprot ID:P12955
Uniprot	P12955
GeneID	5184;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PEPD

Note: This product is for in vitro research use only