

2-oxoglutarate dehydrogenase, mitochondrial

Catalog No: #AP86049



Package Size: #AP86049-1 50ug #AP86049-2 100ug #AP86049-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | 2-oxoglutarate dehydrogenase, mitochondrial |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | 2-oxoglutarate dehydrogenase complex component E1,OGDC-E1,Alpha-ketoglutarate dehydrogenase |
| Accession No. | Uniprot ID:Q02218 |
| Uniprot | Q02218 |
| GeneID | 4967; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OGDH

Note: This product is for in vitro research use only