

MOB kinase activator 1A

Catalog No: #AP86054



Package Size: #AP86054-1 50ug #AP86054-2 100ug #AP86054-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MOB kinase activator 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Mob1 alpha, Mob1A, Mob1 homolog 1B, Mps one binder kinase activator-like 1B, C2orf6, MOB4B, MOBK1B, MOBKL1B
Accession No.	Uniprot ID:Q9H8S9
Uniprot	Q9H8S9
GeneID	55233;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MOB1A

Note: This product is for in vitro research use only