

Neurotrimin

Catalog No: #AP86055



Package Size: #AP86055-1 50ug #AP86055-2 100ug #AP86055-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neurotrimin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hNT,IgLON family member 2,IGLON2,NT,UNQ297/PRO337
Accession No.	Uniprot ID:Q9P121
Uniprot	Q9P121
GeneID	50863;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NTM

Note: This product is for in vitro research use only