

## NAD(P)H-hydrate epimerase

Catalog No: #AP86062



Package Size: #AP86062-1 50ug #AP86062-2 100ug #AP86062-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NAD(P)H-hydrate epimerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Apolipoprotein A-I-binding protein,AI-BP,NAD(P)HX epimerase,Aibp,Apoa1bp
Accession No.	Uniprot ID:B0BNM1
Uniprot	B0BNM1
GeneID	295229;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Naxe

Note: This product is for in vitro research use only