

Methylated-DNA--protein-cysteine methyltransferase

Catalog No: #AP86064



Package Size: #AP86064-1 50ug #AP86064-2 100ug #AP86064-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Methylated-DNA--protein-cysteine methyltransferase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | 6-O-methylguanine-DNA methyltransferase,MGMT,O-6-methylguanine-DNA-alkyltransferase |
| Accession No. | Uniprot ID:Q6TDU1 |
| Uniprot | Q6TDU1 |
| GeneID | 442978; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MGMT

Note: This product is for in vitro research use only