

## DNA mismatch repair protein Msh2

Catalog No: #AP86065



Package Size: #AP86065-1 50ug #AP86065-2 100ug #AP86065-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNA mismatch repair protein Msh2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MutS protein homolog 2
Accession No.	Uniprot ID:P43247
Uniprot	P43247
GeneID	17685;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Msh2

Note: This product is for in vitro research use only