

Mitochondrial peptide methionine sulfoxide reductase

Catalog No: #AP86066



Package Size: #AP86066-1 50ug #AP86066-2 100ug #AP86066-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitochondrial peptide methionine sulfoxide reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Peptide-methionine (S)-S-oxide reductase,Peptide Met(O) reductase,Protein-methionine-S-oxide reductase,PMSR
Accession No.	Uniprot ID:Q9D6Y7
Uniprot	Q9D6Y7
GeneID	110265;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Msra

Note: This product is for in vitro research use only