

Oxidized low-density lipoprotein receptor 1

Catalog No: #AP86069



Package Size: #AP86069-1 50ug #AP86069-2 100ug #AP86069-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oxidized low-density lipoprotein receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Ox-LDL receptor 1, Lectin-like oxidized LDL receptor 1, LOX-1, Lectin-type oxidized LDL receptor 1, Lox1, Oldlr1
Accession No.	Uniprot ID:O70156
Uniprot	O70156
GeneID	140914;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Olr1

Note: This product is for in vitro research use only