

Mannose-binding protein C

Catalog No: #AP86073



Package Size: #AP86073-1 50ug #AP86073-2 100ug #AP86073-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mannose-binding protein C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MBP-C,Collectin-1,MBP1,Mannan-binding protein,Mannose-binding lectin,COLEC1,MBL
Accession No.	Uniprot ID:P11226
Uniprot	P11226
GeneID	4153;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MBL2

Note: This product is for in vitro research use only