

Omega-3 fatty acid receptor 1

Catalog No: #AP86077



Package Size: #AP86077-1 50ug #AP86077-2 100ug #AP86077-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Omega-3 fatty acid receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	G-protein coupled receptor 120,G-protein coupled receptor GT01,Gpr120
Accession No.	Uniprot ID:Q7TMA4
Uniprot	Q7TMA4
GeneID	107221;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:O3far1

Note: This product is for in vitro research use only