

Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial

Catalog No: #AP86079

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86079-1 50ug #AP86079-2 100ug #AP86079-3 1mg

Description

Product Name	Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PDHE1-A type I,PHE1A
Accession No.	Uniprot ID:P08559
Uniprot	P08559
GeneID	5160;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDHA1

Note: This product is for in vitro research use only