

Melanophilin

Catalog No: #AP86081



Package Size: #AP86081-1 50ug #AP86081-2 100ug #AP86081-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanophilin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Exophilin-3,Leaden protein,Slp homolog lacking C2 domains a,Slac2-a,Synaptotagmin-like protein 2a,Ln,Slac2a
Accession No.	Uniprot ID:Q91V27
Uniprot	Q91V27
GeneID	171531;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MIph

Note: This product is for in vitro research use only