

## NADH-ubiquinone oxidoreductase chain 5

Catalog No: #AP86087



Package Size: #AP86087-1 50ug #AP86087-2 100ug #AP86087-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NADH-ubiquinone oxidoreductase chain 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	NADH dehydrogenase subunit 5,MTND5,NADH5,ND5
Accession No.	Uniprot ID:P18940
Uniprot	P18940
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MT-ND5

Note: This product is for in vitro research use only