

NADH dehydrogenase [ubiquinone] iron-sulfur protein 8, mitochondrial

Catalog No: #AP86090

Orders: order@signalwayantibody.com

Package Size: #AP86090-1 50ug #AP86090-2 100ug #AP86090-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	NADH dehydrogenase [ubiquinone] iron-sulfur protein 8, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Complex I-23kD, CI-23kD, NADH-ubiquinone oxidoreductase 23 kDa subunit, TYKY subunit
Accession No.	Uniprot ID:P42028
Uniprot	P42028
GeneID	287027;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NDUFS8

Note: This product is for in vitro research use only