

# NADH dehydrogenase [ubiquinone] iron-sulfur protein 8, mitochondrial

Catalog No: #AP86091

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP86091-1 50ug #AP86091-2 100ug #AP86091-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NADH dehydrogenase [ubiquinone] iron-sulfur protein 8, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Complex I-23kD, CI-23kD, NADH-ubiquinone oxidoreductase 23 kDa subunit
Accession No.	Uniprot ID:Q8K3J1
Uniprot	Q8K3J1
GeneID	225887;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Ndufs8

Note: This product is for in vitro research use only