

## Nuclear receptor subfamily 4 group A member 1

Catalog No: #AP86096



Package Size: #AP86096-1 50ug #AP86096-2 100ug #AP86096-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nuclear receptor subfamily 4 group A member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	NUR77,Nerve growth factor-induced protein I-B,NGFI-B,Orphan nuclear receptor HMR,Hmr,Ngfib
Accession No.	Uniprot ID:P22829
Uniprot	P22829
GeneID	79240;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Nr4a1

Note: This product is for in vitro research use only