

Myoferlin

Catalog No: #AP86097



Package Size: #AP86097-1 50ug #AP86097-2 100ug #AP86097-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myoferlin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Fer-1-like protein 3, Fer1l3, Kiaa1207
Accession No.	Uniprot ID:Q69ZN7
Uniprot	Q69ZN7
GeneID	226101;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Myof

Note: This product is for in vitro research use only