

Dehydrolichyl diphosphate synthase complex subunit Nus1

Catalog No: #AP86101

Orders: order@signalwayantibody.com

Package Size: #AP86101-1 50ug #AP86101-2 100ug #AP86101-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Dehydrolichyl diphosphate synthase complex subunit Nus1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Di-trans,poly-cis-decaprenylcistransferase,Nogo-B receptor,NgBR,Nuclear undecaprenyl pyrophosphate synthase 1 homolog,D10Ert438e,Ngbr
Accession No.	Uniprot ID:Q99LJ8
Uniprot	Q99LJ8
GeneID	52014;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nus1

Note: This product is for in vitro research use only