

NADPH oxidase 3

Catalog No: #AP86108



Package Size: #AP86108-1 50ug #AP86108-2 100ug #AP86108-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	NADPH oxidase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q672J9
Uniprot	Q672J9
GeneID	224480;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nox3

Note: This product is for in vitro research use only