

Protein NKG7

Catalog No: #AP86114



Package Size: #AP86114-1 50ug #AP86114-2 100ug #AP86114-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein NKG7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	G-CSF-induced gene 1 protein,GIG-1 protein,Granule membrane protein of 17 kDa,GMP-17,Natural killer cell protein 7,p15-TIA-1,GIG1
Accession No.	Uniprot ID:Q16617
Uniprot	Q16617
GeneID	4818;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NKG7

Note: This product is for in vitro research use only