

Myosin light chain 3

Catalog No: #AP86115



Package Size: #AP86115-1 50ug #AP86115-2 100ug #AP86115-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Myosin light chain 3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | Myosin light chain 1, slow-twitch muscle B/ventricular isoform,MLC1SB |
| Accession No. | Uniprot ID:P85100 |
| Uniprot | P85100 |
| GeneID | 618352; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYL3

Note: This product is for in vitro research use only