

Neuropeptide FF receptor 2

Catalog No: #AP86116



Package Size: #AP86116-1 50ug #AP86116-2 100ug #AP86116-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuropeptide FF receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	G-protein coupled receptor 74,G-protein coupled receptor HLWAR77,Neuropeptide G-protein coupled receptor,GPR74,NPFF2,NPGPR
Accession No.	Uniprot ID:Q9Y5X5
Uniprot	Q9Y5X5
GeneID	10886;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NPFFR2

Note: This product is for in vitro research use only