

## Nuclear receptor subfamily 4 group A member 2

Catalog No: #AP86118



Package Size: #AP86118-1 50ug #AP86118-2 100ug #AP86118-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nuclear receptor subfamily 4 group A member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	NUR-related factor 1,Nuclear orphan receptor HZF-3,Orphan nuclear receptor NURR1,Regenerating liver nuclear receptor 1,RNR-1,SL-322,Hzf-3,Nurr1,Rnr1
Accession No.	Uniprot ID:Q07917
Uniprot	Q07917
GeneID	54278;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Nr4a2

Note: This product is for in vitro research use only