

Metaxin-1

Catalog No: #AP86120



Package Size: #AP86120-1 50ug #AP86120-2 100ug #AP86120-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Metaxin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mitochondrial outer membrane import complex protein 1,Mtx,Mtxn
Accession No.	Uniprot ID:P47802
Uniprot	P47802
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mtx1

Note: This product is for in vitro research use only