

Proton myo-inositol cotransporter

Catalog No: #AP86121



Package Size: #AP86121-1 50ug #AP86121-2 100ug #AP86121-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proton myo-inositol cotransporter
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	H(+)-myo-inositol cotransporter,H(+)-myo-inositol symporter,Solute carrier family 2 member 13
Accession No.	Uniprot ID:Q3UHK1
Uniprot	Q3UHK1
GeneID	239606;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc2a13

Note: This product is for in vitro research use only