

2'-5'-oligoadenylate synthase 3

Catalog No: #AP86124

Package Size: #AP86124-1 50ug #AP86124-2 100ug #AP86124-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2'-5'-oligoadenylate synthase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	(2-5')oligo(A) synthase 3
Accession No.	Uniprot ID:Q5MYT7
Uniprot	Q5MYT7
GeneID	494202;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Oas3

Note: This product is for in vitro research use only