

Protein kinase C-binding protein NELL2

Catalog No: #AP86125



Package Size: #AP86125-1 50ug #AP86125-2 100ug #AP86125-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein kinase C-binding protein NELL2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	NEL-like protein 2, Nel
Accession No.	Uniprot ID:Q62918
Uniprot	Q62918
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nell2

Note: This product is for in vitro research use only