

Neuromodulin

Catalog No: #AP86126



Package Size: #AP86126-1 50ug #AP86126-2 100ug #AP86126-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuromodulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Axonal membrane protein GAP-43, Calmodulin-binding protein P-57, Growth-associated protein 43, Basp2
Accession No.	Uniprot ID:P06837
Uniprot	P06837
GeneID	14432;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gap43

Note: This product is for in vitro research use only