

Nuclear factor NF-kappa-B p105 subunit

Catalog No: #AP86127



Package Size: #AP86127-1 50ug #AP86127-2 100ug #AP86127-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear factor NF-kappa-B p105 subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
Accession No.	Uniprot ID:Q6F3J0
Uniprot	Q6F3J0
GeneID	442859;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NFKB1

Note: This product is for in vitro research use only