

## Substance-P receptor

Catalog No: #AP86129



Package Size: #AP86129-1 50ug #AP86129-2 100ug #AP86129-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Substance-P receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SPR,NK-1 receptor,NK-1R,Tachykinin receptor 1,Tac1r
Accession No.	Uniprot ID:P30548
Uniprot	P30548
GeneID	21336;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tacr1

Note: This product is for in vitro research use only