

Diphosphoinositol polyphosphate phosphohydrolase 3-beta

Catalog No: #AP86136

Orders: order@signalwayantibody.com

Package Size: #AP86136-1 50ug #AP86136-2 100ug #AP86136-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Diphosphoinositol polyphosphate phosphohydrolase 3-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DIPP-3-beta,Diadenosine 5',5'''-P1,P6-hexaphosphate hydrolase 3-beta,Diadenosine hexaphosphate hydrolase (AMP-forming),Nucleoside diphosphate-linked moiety X motif 11,Nudix motif 11,Dipp3b,MNCb-1696
Accession No.	Uniprot ID:P0C028
Uniprot	P0C028
GeneID	102954;58242;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nudt11

Note: This product is for in vitro research use only