

## Osteopontin

Catalog No: #AP86137



Package Size: #AP86137-1 50ug #AP86137-2 100ug #AP86137-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Osteopontin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Bone sialoprotein 1, Secreted phosphoprotein 1, SPP-1, 2b7, Spp-1
Accession No.	Uniprot ID:P08721
Uniprot	P08721
GeneID	25353;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Spp1

Note: This product is for in vitro research use only