

Ubiquitin thioesterase otulin

Catalog No: #AP86138

Package Size: #AP86138-1 50ug #AP86138-2 100ug #AP86138-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquitin thioesterase otulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Deubiquitinating enzyme otulin, OTU domain-containing deubiquitinase with linear linkage specificity, Ubiquitin thioesterase Gumby, Fam105b, Gum
Accession No.	Uniprot ID:Q3UCV8
Uniprot	Q3UCV8
GeneID	432940;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Otulin

Note: This product is for in vitro research use only