

## Oxysterols receptor LXR-alpha

Catalog No: #AP86143



Package Size: #AP86143-1 50ug #AP86143-2 100ug #AP86143-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Oxysterols receptor LXR-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Liver X receptor alpha,Nuclear receptor subfamily 1 group H member 3,LXRA
Accession No.	Uniprot ID:Q13133
Uniprot	Q13133
GeneID	10062;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NR1H3

Note: This product is for in vitro research use only