

Protein orai-3

Catalog No: #AP86148



Package Size: #AP86148-1 50ug #AP86148-2 100ug #AP86148-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein orai-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Transmembrane protein 142C,Tmem142c
Accession No.	Uniprot ID:Q6AXR8
Uniprot	Q6AXR8
GeneID	309000;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Orai3

Note: This product is for in vitro research use only