

A-Myb Antibody

Catalog No: #44287

Package Size: #44287-1 50ul #44287-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	A-Myb Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	A-Myb Antibody detects endogenous levels of total A-Myb
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human A-Myb
Target Name	A-Myb
Other Names	MYBL1, A-MYB, Myb-like protein 1, AMYB, Myb-related protein A
Accession No.	Swiss-Prot#: P10243NCBI Gene ID:4603
Uniprot	P10243
GeneID	4603;
Calculated MW	85kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of A-MYB expression in mouse brain (A), rat brain (B) whole cell lysates.

Note: This product is for in vitro research use only