

N-acetylaspartate synthetase

Catalog No: #AP86158



Package Size: #AP86158-1 50ug #AP86158-2 100ug #AP86158-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-acetylaspartate synthetase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NAA synthetase,Camello-like protein 3,N-acetyltransferase 8-like protein,CML3
Accession No.	Uniprot ID:Q8N9F0
Uniprot	Q8N9F0
GeneID	339983;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NAT8L

Note: This product is for in vitro research use only