

Myosin light chain 3

Catalog No: #AP86162



Package Size: #AP86162-1 50ug #AP86162-2 100ug #AP86162-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin light chain 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Myosin alkali light chain 1, ventricular/slow skeletal muscle isoform, Myosin light chain 1, slow-twitch muscle B/ventricular isoform, MLC1SB, Mlc1v, Mylc
Accession No.	Uniprot ID:P09542
Uniprot	P09542
GeneID	17897;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Myl3

Note: This product is for in vitro research use only