

## Homeobox protein NANOG

Catalog No: #AP86163



Package Size: #AP86163-1 50ug #AP86163-2 100ug #AP86163-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Homeobox protein NANOG
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ES cell-associated protein 4, Early embryo specific expression NK-type homeobox protein, Homeobox transcription factor Nanog, Ecat4, Enk
Accession No.	Uniprot ID:Q80Z64
Uniprot	Q80Z64
GeneID	71950;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Nanog

Note: This product is for in vitro research use only