

Mimecan

Catalog No: #AP86166



Package Size: #AP86166-1 50ug #AP86166-2 100ug #AP86166-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mimecan
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Osteoglycin,OIF
Accession No.	Uniprot ID:P19879
Uniprot	P19879
GeneID	280884;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OGN

Note: This product is for in vitro research use only