

## Poly [ADP-ribose] polymerase 2

Catalog No: #AP86168



Package Size: #AP86168-1 50ug #AP86168-2 100ug #AP86168-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Poly [ADP-ribose] polymerase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PARP-2,ADP-ribosyltransferase diphtheria toxin-like 2,ARTD2,NAD(+) ADP-ribosyltransferase 2,ADPRT-2,Poly[ADP-ribose] synthase 2,pADPRT-2,Adprt2,Adprt12,Aspart12
Accession No.	Uniprot ID:O88554
Uniprot	O88554
GeneID	11546;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Parp2

Note: This product is for in vitro research use only